


# 1600 MiniG™

Mixing and  
Blending 

Grinding and  
Pulverizing 

## POWERFUL, MINI, TISSUE HOMOGENIZER



The MiniG™ is the ideal solution for labs that want the power of the 2010 Geno/Grinder® with a lower sample throughput. It was created for applications involving sample preparation for DNA/RNA and protein extractions as well as pesticide residue analysis.

### TYPICAL SAMPLES

Plant and animal tissue, cell cultures, seeds, yeast, leaves and bacteria.

### FEATURES

- • Vigorous clamp movement results in increased yields.
- • Consistent vertical motion of the clamp ensures equal conditions for all samples.
- • Kryo-Tech® accessories are available to preserve temperature sensitive samples.
- • Light & portable with a small footprint that minimizes bench space.

SPEX® SamplePrep® 

[www.spexsampleprep.com](http://www.spexsampleprep.com)



# 1600 MiniG™

## Specifications

1. Window in lid allows you to see the grinding process.
2. Lid has a safety interlock switch that stops the motor if the lid is opened during the grinding process.
3. Easy-to-use clamp keeps samples secure during the shaking process.
4. Clamp typically holds up to 2 titer plates or 48 x 2 mL vials when stacked.
5. LCD display screen shows run times in minute: second format.
6. Rate and timer controls allow you to adjust speed and set time.
7. Start and stop buttons allow you to start and stop the grinding process.

<b>Voltage</b>	115V/60Hz or 230V/50Hz	<b>Power Cord</b>	3-prong grounded plug 115V / 60Hz. 2-prong European cord for 230V / 50Hz.
<b>CE Approved</b>	Yes		
<b>Dimensions</b>	21.5 in. (54.6 cm) high x 11.25 in. (28.6 cm) wide x 17.0 in. (43.2 cm) deep	<b>Clamp Movement</b>	1.25 in. (3.2 cm) vertical
<b>Weight</b>	40 lbs. (22 kg)	<b>Clamp Speed</b>	Adjustable from 500 to 1500 strokes/minute
<b>Motor</b>	1/7 HP	<b>Timer/Controls</b>	Digital display in minutes & seconds, maximum 10:00

### Corporate Headquarters

SPEX SamplePrep  
15 Liberty Street  
Metuchen, NJ 08840  
1-855-GET-SPEX

Tel: 732-623-0465  
Fax: 732-906-2492  
sampleprep@spex.com  
www.spexsampleprep.com

### European Division

SPEX CertiPrep Ltd  
2 Dalston Gardens  
Stanmore, Middlesex HA7 1BQ

Tel: +44(0)208-204-6656  
Fax: +44(0)208-204-6654  
sales@spexcertiprep.co.uk  
www.spexcertiprep.com