

DigiPREP DISTILLATION SYSTEM

★ Systems ranging from manual to fully automatic

- Fast steam distillation system
 - 3 - 5 minutes
- Automatic shutoff
- pH titration (optional)


★ Applications

- Kjeldahl, TKN, Ammonia
- Alcohol in wine and beer
- Phenol, Cyanide & Fluoride



FEATURES

	100	200	300	400	500
Automatic addition of H ₂ O			•	•	•
Automatic addition of NaOH		•	•	•	•
Automatic addition of H ₃ BO ₃				•	•
Programmable reaction time		•	•	•	•
Programmable distillation time		•	•	•	•
Automatic steam generation	•	•	•	•	•
Variable steam generating capacity (40% -100 %)		•	•	•	•
Automatic removal of sample waste by suction			•	•	•
Number of programs	0	1	1	99	99
Display language can be chosen by the user		•	•	•	•
Visual Error Reports	•	•	•	•	•
Acoustic Error Reports		•	•	•	•
Serial Interface (RS232)			•	•	•
Standby operation between distillations	•	•	•	•	•
Level monitoring for the set of canisters	•	•	•	•	•
Connection for titration unit					•



SCP SCIENCE

www.scpscience.com

Sample
Prep

ICP-AES & MS

AA

XRF

Calibration
Standards

Certified
Reference
Materials

Also Available

*DigiPREP HT
High Temperature
Digestion System*



The *DigiPREP HT* is a new System for Kjeldahl digestion
 -Versatile heating to 450°C
 -Acid resistant coated graphite block
 -Available in two sizes:
 40 tube (100ml) or
 20 tube (250ml) sample capacity

Call for more information or to request a quotation!!

MK-MKG003-DIST-1.0-E

USA

Tel.: (800) 361-6820
 Fax: (800) 253-5549

Canada / International

Tel.: (800) 361-6820 / (514) 457-0701
 Fax: (800) 253-5549 / (514) 457-4499

Europe

Tel.: +33 (0)1 69 18 71 17
 Fax: +33 (0)1 60 92 05 67